

Solid Zirconia Crowns - Solid Dental Zirconia Crowns for bruxers - LVDDS Dental Zirconia Lab


[HOME](#) [COMPANY](#) [SHIPPING](#) [CONTACT US](#)
[BITE SPLINTS](#) [FULL-CAST & PFMS](#) [IMPLANT ABUTMENTS](#) [METAL-FREE](#)

BruxZir® Solid Zirconia

Solid zirconia restoration with no porcelain overlay
 Virtually unbreakable zirconia is ideal for bruxers and grinders
 Prescribe instead of metal PFMs and cast gold
 Smooth surface wears gently on opposing teeth

[BRUXZIR SOLID ZIRCONIA](#) [IPS E.MAX CAD](#) [PRISMATIK CZ](#) [3M ESPE LAVA](#)

BruxZir® Solid Zirconia

DENTISTS

BruxZir® is a solid zirconia crown or bridge dental restoration with no porcelain overlay. More brawn than beauty, you will be impressed by the esthetics of BruxZir when prescribed instead of metal occlusal PFM and full-cast gold dental restorations. Virtually chip-proof, BruxZir solid zirconia crowns and bridges are the ideal dental restoration for your bruxing and grinding patients.

Designed and milled using advanced CAD/CAM technology, BruxZir zirconia crown and bridge restorations are sintered for 6.5 hours at 1,530 degrees Celsius. The final BruxZir solid zirconia crown or bridge emerges nearly "bulletproof" and is glazed to a smooth surface.

To send a case to our dental lab: Conservatively prepare as thin as 0.5 mm with feather edge margins, much like you would cast gold. Any preparation with a minimum 0.5 mm of occlusal space is acceptable, but 1 mm is ideal.

To learn more about the benefits of BruxZir, visit www.bruxzir.com.

LAB OUTSOURCE

BruxZir is a solid zirconia crown or bridge restoration with no porcelain overlay. More brawn than beauty, you will be impressed by the esthetics of BruxZir when prescribed instead of metal occlusal PFM and full-cast gold crowns restorations. Virtually chip-proof, BruxZir solid zirconia crowns and bridges are the ideal dental restoration for your bruxing and grinding patients.

Designed and milled using CAD/CAM technology, BruxZir zirconia crown and bridge restorations are sintered for 12 hours at 1,530 degrees Celsius. The final BruxZir solid crown or bridge emerges nearly "bulletproof" and is glazed to a smooth surface.

To send a case to our dental lab: Send modelwork or 3Shape digital design file. Indicate final ceramic shade on lab work order.

PRODUCT IMAGES



[Metal-Free Preparation Guide](#)



[Rx for Dentists](#)



[Rx for Lab Outsource](#)

[HOME](#) [COMPANY](#) [SHIPPING](#) [HIPAA](#) [SITE MAP](#) [CONTACT US](#) [Taxpayer ID #95-2743032 \(W-9 PDF\)](#)

For questions or assistance, call us at 800-936-1848 (M-F 7:00AM - 5:00PM PST)
 © 2012 Las Vegas Digital Dental Solutions