

Exhibit 89

Assured Dental Lab

Quality Products at Amazing Prices.

Case after case.

Home

About Us

Products & Services

Customer Service

News, Events & FAQs

Contact Us

Resources/Manufacturers

Site Map

P.O. Box 4420 Portland, OR 97208-4420

Toll free: 1-877-283-5351

Products & Services

Download our [RX Forms](#) to get started today or call us at **1-877-283-5351** for additional product information.

Dental Product Information



Full implant solutions for all indications.

We work with a specialized team, with years of experience in the design and manufacturing of dental implants. We offer diverse and the most popular [implant systems](#).



Duraflex™



Duraflex Partial Denture delivers patient satisfaction!

- Natural looking-virtually disappears in the mouth
- Easy to adjust, even easier to polish and it's repairable
- Won't absorb moisture: stain and odor resistant
- Reliable-virtually unbreakable
- Comfortable thermoplastic material-lightweight and flexible



Calypso™



The New Generation All-Ceramic Restorations are supported by a high strength zirconia framework and given the translucency nature of zirconia, all restorations closely mimic the appearance of natural dentition. The Calypso Zirconia framework is available in 8 vita classic core shades to optimize color and clarity. Zirconia, at 1200 MPa flexural strength, is better able to resist crack propagation and with two overlays of porcelain, the restorations have life-like aesthetics. Indications are anterior, posterior and up to 6-unit bridges.



FMZ is Now Available!



FMZ crowns and bridges are full-milled, pure Zirconia which means no porcelain and no possibility of porcelain fractures! As a full gold crown replacement, FMZ offers: natural aesthetics and lifelike characteristics, 1200 MPa flexural strength, adjustability and stainability, excellent marginal fit, conventional chamfer prep and cementation. FMZ is available at half the price of a gold crown. Indications are: crowns and up to 4-unit bridges, crowns over implants and crowns with limited occlusal clearance.



3M™ ESPE™ Lava™



3M ESPE Lava Crown and Bridges are a part of the Lava Precision Solutions System for anterior and posterior crowns and bridges up to 48mm in length. Lava™ Crowns and Bridges delivers: precision CAD/CAM system for outstanding marginal fit, proven high-strength performance-1440 MPa, natural translucency and patented shading system available in 8 shades producing unmatched aesthetics. Lava restorations may be cemented with traditional or self-adhesive resin cement and prep requires shoulder or chamfer to allow precise margin detection. Approved indications are: single crowns; 3- and 4-unit bridges with multiple abutments; 5- and 6-unit bridges; Splinted crowns (up to 6 units); Cantilever bridges (up to 4 units); custom implant abutments.



Metal Cast Partial



Metal Cast Partial represent an excellent practical and affordable prosthetic choice available for your patients. Even though cast metal partials are made of stronger materials, they are less bulky. They fit better and last longer, providing years of comfort for your patient. When appropriate, their use should be considered as a viable treatment option. We offer a wide variety of removable prosthetics featuring:

- Non-Allergenic (nickel and beryllium free)
- Unsurpassed memory
- Ease of adjustment
- Beautiful polishability
- Using Ivoclar Vivadent Premium BlueLine Teeth



Dentures

It is clear that there still is a tremendous need for removable appliances today. In the United States alone, it



is estimated that there will be approximately 38 million edentulous elderly adults by the year 2020 compared with only 33.6 million in 1991; the aging of the "baby boom" population largely accounts for the increase. While the basic process of making dentures has changed little over the past several decades, new materials and techniques provide functional, esthetic restorations that offer exceptional value to patients. Implant treatment is a tremendous adjunct to removable prosthodontics in the treatment of edentulous patients, but it is not within the financial reach of all dental patients. Assured Dental Laboratory possesses the knowledge and trained technicians to deliver complete dentures that mimic the color and characterization of natural gums. All dentures come in a variety of shades with Ivoclar Vivadent Premium BlueLine teeth.



Ceramage®

Achieve PURE esthetic results with Ceramage, a Zirconium Silicate integrated indirect restorative with positive 5-year clinical studies. A PFS (progressive fine structure) filling of more than 73% plus an organic polymer matrix delivers superior flexural strength and unsurpassed polishability. Ceramage® provides superior wear characteristics; simulating the wear and translucency of natural dentition, excellent for implants, inlays, onlays, crowns and bridges. Ceramage allows for conventional cementation and is extraordinarily biocompatible.



IPS e.max

IPS e.max comprises highly esthetic and high strength materials for both the press and the cadcam techniques. Recommend for **single units & up to 4-unit bridges** conventional preparation and cementation. With e.max you can choose the framework material. Either glass ceramic or zirconium oxide the new material generation is based on nano-fluorapatite and endows the restorations with utmost individual appearance and true to nature properties.



IPS Empress Esthetic

IPS Empress Esthetic raises the bar for **pressable ceramics**. For over 15 years Empress still sets the standard for all ceramic esthetics. Offering:

- » True to nature esthetics
- » Homogeneous leucite crystal distribution allowing for higher flexural strength
- » Improved optical properties
- » Optimized color and clarity
- » 11 Ingot shades extending shades in bleach range



Captek

The ultimate in **porcelain fused to metal crown and bridge**, Captek offers all ceramic esthetics and fit. Recommend for single crowns and up to 3 unit bridges. The capillary technology offers dentist and patients the beauty, health and longevity of gold produced through a new process with unparallel strength accuracy and esthetic. [Click here for more information.](#)



Cercon Zirconia

Smart Ceramics for Crowns and Bridges

Yttria-stabilized Zirconia oxide is an innovative system for **metal-free single crowns and up to 6-unit bridges**. Recommended for any placement in the mouth, Cercon is the strongest and the toughest **all-ceramic restoration** being offered today. Conventional preparation, cementation and adhesive bonding techniques. Over 2,000,000 units placed since its introduction.

Home :: About Us :: Products & Services :: Customer Service :: News, Events & FAQs ::
Contact Us :: Resources/Manufacturers :: Site Map

Exhibit 90

[Home](#) [About Us](#) [Lab Tour](#) [Services](#) [Restorative Options](#) [Products](#) [Testimonials](#) [Seminars](#) [Print Rx/Labels](#) [Clinic Video/Links](#) [Shipping/Schedule](#) [Contact](#)

Restorative Options

Authentic Dental Laboratory is proud to be a full service dental laboratory. Our Laboratory features some of the most current, technically advanced products that deliver reliable restorations with superior esthetics at competitive prices.

"NEWS FLASH"

Authentic Lab owns and operates the most advanced CAD/CAM milling center in South Texas. Call us for your milling and digital scan needs

LAVA



Lava offers precision CAD/CAM milling for outstanding marginal fit. Zirconia's exceptional strength has 5 years of clinical history with 100% survival rate. Natural translucency and excellent esthetics. A good choice for singles anywhere and bridges up to 12 units.

Excellent Fit
Strength
Esthetics
Proven Performance

Procera combined with All-Ceram porcelain allows us to recreate the beauty of natural teeth. You can cement with the crown & bridge cement of your choice or bond in place. Because of Procera's exceptional strength utilizing aluminium and zirconia cores you can prescribe it for any position in the mouth with confidence that it will not fracture under occlusal loading. CAD-CAM technology makes the passive fit possible. Once again, a 360° 1-1.5mm chamfer prep with rounded internal line angles and 1.5-2mm occlusal reduction is protocol.

Procera



Procera Bridge

Procera Bridge is an all-ceramic restorative system from

Ceramco that delivers the strength of a PFM and the beauty of all-ceramics. Bridges up to 14 units anywhere.



Empress has been a successful and widely-used all-ceramic system. This pressed or milled, highly translucent ceramic has beautiful, life-like esthetics and wears like natural teeth. Empress is indicated for anterior single crowns and veneers and must be adhesively bonded.

Empress



Emax



Emax offers the possibility of all-ceramic bridges to the 2nd bicuspid. Super-strong lithium disilicate frameworks and the fluorapatite layering ceramic allow for restorations very similar to natural teeth. Bond or cement – your choice.

ChromaTech Veneers exclusively from Authentic Dental Lab are fabricated using the most advanced porcelain building techniques to achieve polychromatic esthetics. A top choice for many cosmetic dentists, ChromaTech was developed specially for esthetics in the anterior region of the mouth.

ChromaTech Veneers



GC Gradia



GC Gradia is a unique light-cured Micro-Ceramic-Composite system for crowns, bridges and veneers. The result of this restoration has an opalescence compatible to porcelain and close to natural dentition.

Captek is a 22k gold substructure that enhances the porcelain giving a soft, warm, natural look. Capillary technology makes a better fit possible. The high gold content of Captek makes it an excellent choice for biocompatibility. The finest restorative system available for porcelain fused to metal crowns or bridges. Highly recommended because of esthetics, strength and purity of materials. Captek is truly the leading PFM of choice.

Captek



PFM



PFM patients today do not want pretty crowns, they want beautiful teeth. Utilizing the most advanced porcelains today that contain opalescence and natural fluorescence and applied by a skilled ceramist makes crowns look like teeth.

Implants at Authentic Dental Laboratory are supported by the team approach to implant dentistry. As a result, we can provide you with a diagnostic case plan to assist you in case design and planning, pre and post surgical treatment and a cost estimation on the proposed reconstruction.

Implants



Partials / Dentures / Attachments

Partials / Dentures / Attachments Authentic Lab's removable department offers the best in quality, service, technology, timeliness and esthetics at the prices you and your patients will appreciate. Premium



partials, dentures, comfort night guards, flippers, attachment cases and more with superior personalized service.

Milled Bars are an excellent option today when the placement of implants are a challenge. These bars are milled to perfection, and as you can see, there are no limits even when restoring more difficult cases.

Milled Bars



Snap On



Snap-On smiles are an inexpensive alternative to conventional procedures. Great for implant temporaries, cosmetic enhancement and to replace missing teeth. Most patients use this appliance as a stepping stone to permanent restorative dentistry.

Zirconia Full Contour

Zirconia Full Contour Crowns are available in either Lava or Authentic Zirconia. This monolithic crown is great for bruxing patients, replacing full metal crowns or metal occlusion in posterior bridges or single units. With a (Lifetime Warranty) you can't go wrong.



Hybrid Zirconia




Hybrid Zirconia crowns and bridges are an excellent product when occlusal clearance is a challenge. Zirconia is placed on the occlusal surface then veneered with porcelain to achieve excellent esthetics.



1950 Bandera Rd. San Antonio, Texas 78228
800-683-1025 210-735-1433 FAX: 210-735-2127
services@authenticlab.com

Exhibit 91



PRODUCTS

CHINA DENTAL OUTSOURCING

Quality is our motto

PRODUCTS

All Ceramic: All Zirconia for Bruxers

Prosthetics - All Zirconia for Bruxers

All Zirconia for Bruxers is an innovative system for metal-free crowns and bridges in the posterior used primarily for bruxers and heavy biters. The strength of zirconia restorations is well established, and the ability to provide an all zirconia restoration. The system is used to fabricate anterior or posterior single units or multi-unit bridge restorations out of translucent, biocompatible zirconia, the strongest and toughest dental ceramic.

Esthetically, they are inferior to porcelain restorations, but are superior to full metal or metal occlusal. Shades are close approximations (but cannot be exact) to the Vita 16 shade guide, hence they are best used in the posterior.

They are extremely strong; being crack and fracture resistant, and extremely durable.

[Contact us for a Price Quote.](#)

PRODUCTS

[Crown and Bridge](#)
[Captek](#)
[All Ceramic - E.max](#)
[All Zirconia for Bruxers](#)
[Circon Zirconia](#)
[Veneers & e.Max Thin Veneers](#)
[Composites](#)
[Full Dentures](#)
[Partial Denture Framework](#)
[Valplast Flexible Dentures](#)
[Implants](#)
[Other Types](#)
[Lab Photos](#)
[Brands](#)
[List of Lab Services](#)

Enlarged Photos to Show Detail



FOR DENTAL LABS & DENTISTS ONLY
Have Questions?
Call us in China!



USA Toll-free
 1 888-262-8583

Australia Toll-Free
 1800-557-863

[International Phone Codes for All Countries](#)

Current Time in China is:
Tue, Nov 27, 2012
3:04:27 AM

We answer in China from 9 AM to 11 PM China time 7 days a week!

We answer in the U.S. from 9 AM to 5:30 PM CST time Monday-Friday

Or call anytime and leave a message!

[FOLLOW ME ON twitter](#)



You and 69 others like this. 69 people like this. Be the first of







[Contact us for a Price Quote.](#)



[HOME](#) [PRODUCTS](#) [TESTIMONIALS](#) [DELIVERY TIME](#) [PRICE QUOTE](#) [POLICIES](#) [FAQs](#) [ABOUT US](#) [CONTACT US](#)

Copyright © China Dental Outsourcing, Inc., 1998-2012. All Rights Reserved


Exhibit 92

- [Porcelain Fused to Metal](#)
- [Full Cast Restorations](#)
- [CAD/CAM & Non Metallics](#)
 - [Lava](#)
 - [Procera](#)
 - [IPS e.max](#)
 - [Cerec Connect](#)
 - [BruxZir](#)
- [Removables](#)
- [Composites & Temporaries](#)
- [Implants](#)

★ **How to Get Started**

1. [Request Intro Kit](#)
2. [Print RX](#)
3. [Schedule Case Pickup](#)

Lava Prep & Handling Information



Lava Prep Guide

[Download File \(2MB\)](#)



CAD/CAM & Non-Metallic
Authorized Lava Milling Center



Specifications:

- 1440 MPa Proven High Strength
- 5 year clinical study with 100% survival rate
- Patented substructure available in 8 different shades
- Requested 3X more than the leading brand
- Comes with barcoded authenticity label

Indications:

Single unit crowns and bridges


Cementation:

For conventional cementation - use a glass ionomer cement (i.e. Ketec™ Cem) or for self adhesive cementation please use 3M RelyX Unicem.

For additional information, please see page 14 in the Prep and Handling Guide (on the left)

Our dental laboratory specializes in **lava crowns** and **lava bridges**. We are an authorized **3M ESPE Lava Milling Center**. That means each lava crown is barcoded for authenticity and is made from the finest zirconia available. **Lava Zirconia** restorations are an ideal alternative to porcelain-fused-to metal (PFM) restorations. The zirconium oxide achieves the beauty and translucency of natural teeth. Over time, metal restorations can leave a dark metal margin where the gums recede. With **Lava zirconia**, it produces strong durable restorations and retains its original beauty for years. It is a much better alternative for those patients who have concerns about metal.

Lava™ is a registered trademark of 3M ESPE

**Customer Poll**

Which all-ceramic restoration do you prefer?

☐ IPS E.max CAD

☐ 3M ESPE Lava

☐ Procera

☐ IPS Empress

☐ Generic zirconia

☐ Other

Vote

Results





fig.2 Light intensity transmitted through zirconia copings as detected by photo-radiometer

▲ Go to top

Lava Authenticity



Now you can verify the authenticity of your lava restorations. Just go to 3M's website and enter your unique lava barcode.

Lava Chairside Oral Scanner



As an authorized Lava C.O.S. laboratory, you can offer your patients the convenience and comfort of digital impressions

ST Design Center



Lava design centers can enjoy flat-rate pricing and convenient 24-48hr turnaround time for their zirconia copings.



Exhibit 93



Matt 11:28-30

- Products
- Pickup and Delivery
- Dentist Support & Forms
- My Account
- Contact Us
- Careers
- About PDL
- Home



Downloads

- [Print Standard Rx](#)
- [Print Proline Rx](#)
- [Request Supplies](#)

Contact

Pittman Dental Laboratory (PDL)
2355 Centennial Circle
Gainesville, GA 30504-5799 USA
Local: (770) 534-4457
Free: (800) 235-4720
Fax: (770) 503-1173
[Click to E-Mail](#)
www.pittmandental.com

Turn Around Times in Working Days (including shipping)	
PFM(up to 4 units)	10
Full Cast	10
All-Ceramic	12
Captek	12
Bridges(5 units & up)	12

Dentures and Partials

Framework	11
Setup/Tryin	8
Finish	8

[Request more information.](#)

Call Toll Free: +(800) 235-4720

ZIRCONIA

Pittman Dental is proud to offer the following All-Ceramic Crown:

- **Lava**
- **Procera Zirconia**

Advantages/Patient Benefits

- High Strength And Durability
- Bio-Compatible With High Aesthetics
- Only Reliable Alternative For Patients With Metal Allergies

Cementation

- Traditional Cementation Or Adhesive Bonding May Be Used.

Indications

- Crowns
- Bridges: Of All Kinds And Different Lengths
- Implants: Abutments And Crowns

Prep Guide

Call 800-235-4720 for information.

Made by Express Marketing Solutions LLC,

Copyright © Pittman Dental Laboratory, Inc. (PDL)
All Rights Reserved: all trademarks and brands are property of their respective owners.
[Privacy Policy](#)

Exhibit 94



Full Contour Zirconia Crowns



Our **Bruxzir and Pearl crowns** are a durable solution that is less expensive and more esthetically pleasing than full cast crowns. Because frequent grinding or clenching of teeth puts extreme wear on dental work, your bruxing patients need a strong solution. Our Bruxzirs and Pearls deliver, with an all-zirconia structure that provides material durability... without the cost of typical PFM restoration or the obvious look of full cast metal crowns. The CAD/CAM milling technology delivers consistent quality while the all-zirconia material provides solid, chip-resistant results with no porcelain overlay.

INSTRUCTIONS FOR ADJUSTMENT AND POLISHING

Adjust using water and air spray to keep the restoration cool and to avoid micro-fractures with a fine grit diamond. If using air only, use the lightest touch possible when making adjustments. A football shaped bur is most effective for occlusal and lingual surfaces on anterior teeth; a tapered bur is the ideal choice for buccal and lingual surfaces.

Somer Dental recommends using a Brasseler USA Zirconia adjustment kit. To order call Brasseler @ 1-800-841-4522 : product order # 225.00.638. For a modified prep kit order # K000017.

INSTRUCTIONS FOR CEMENTATION

- o Resin-reinforced glass ionomer cement (Rely X Luting Cement, 3M ESPS, GC Fuji Plus, GC America)
- o Resin cements for short or over-tapered preparations (Rely X Unicem, 3M ESPE)

This post was written by:

Elvis - who has written 18 posts on Somer Dental Labs.

[Contact the author](#)

Leave a Reply

 Name (required)

 Mail (will not be published) (required)

 Website


Somer Dental Laboratory was founded in 1982 by Larry Sowinski, CDT and Keith Spencer, CDT. Somer is located in Zionsville, Indiana.

Everyday, our technicians set teeth, build porcelain, cast frames, wax copings, scan models, and mill crowns on schedule on time and get great report cards from our dentists. [Learn More](#)



certified
dental
laboratory

We are proud to be a Certified Dental Lab. Certification enables us to show you we adhere to the highest standards of quality control, materials and customer service.



Rx Form

[Click here to download our RX Form.](#)

Somer Hours of Business

We are open for business and welcome your inquiries from 7:00 a.m. to 5:00 p.m. EDT, Monday through Thursday. Our drivers pick-up and deliver cases Monday through Friday.

To schedule a pick-up, call 317.873.1111 or 800.283.1110 during business hours or Friday mornings. After hours, please call the Somer Pick-up Hotline at 317.873.1119.

Or [Email Us](#) for more information.

Somer's Vision and Mission Statement

Vision:

In unison, as a Team, Somer Dental Laboratories' vision is to establish an environment, which can productively provide quality products and services, education, personal growth, and financial security.

Submit Comment

« Previous Entries

Next Entries »

Mission Statement:

Somer Dental Laboratories is committed to be a leader in the dental laboratory industry by educating our team to provide the highest standards of dental care through excellent service, quality products and advanced technology.

DDX Dental

Sign Up to submit online cases with Somer Dental Laboratories & DDX

Go to Somer Dental Laboratories's DDX dashboard

Go to Somer Dental Laboratories DDX Rx



Exhibit 95

TRACHSEL DENTAL STUDIO

HomeAbout UsProductsResourcesSend a CaseContact

<

NADL
National Association of
Dental Laboratories

CERTIFIED
DENTAL
LABORATORY

>

About UsProductsContact Us

Trachsel Dental Studio is a full service dental lab located in Rochester, MN. TDS is a proud member of the National Dental Lab Association (NADL) and the Minnesota Dental Laboratory Association (MDLA).

Trachsel Dental has 35 employees including four full time Certified Dental Technicians. Give us a call or send us your next case, we guarantee you'll be impressed!

[Learn more about us](#)

Crown & Bridge

BruxZir®

Captek™

Porcelain Fused to Metal

Fortress™

Pressed on Metal

Gold Inlays/Onlays

Porcelain to Zirconia

Full Cast Restorations

IPS e.max® Press

Dentures & Partial

Orthodontics

Implants

name

email

message

SEND

NavigationLatest TweetContact Us

Home

About Us

Products

Resources

Send a Case

Contact

Have you experienced full contour zirconia..? Were an authorized BruxZir lab, guaranteed youll be impressed! <http://twitpic.com/52f3ib>

t

Follow Us

f

Find Us

Trachsel Dental Studio
1834 15th St. NW
Rochester, MN 55901

(Toll) (800) 831-2362
(Ph) (507) 288-2362
(F) (507) 288-3972

© 2011-2012 Trachsel Dental Studio | [Zennic](#) | [Login](#)

Exhibit 96



Fixed

[IPS e.max® Press](#)

[Procera AllCeram](#)

[Lava™](#)

[IPS Empress®](#)

[Sculpture®/Fibrekor®](#)

[BruxZir](#)

[Captek](#)

[BIO 2000®](#)

[Conventional PFMs](#)

Removables

Implants

Radiographic & Surgical

Appliances





York Dental / Cerama Laboratory, Inc. © 2010

Designed by AMG Creative

Exhibit 97



SHOWCASE
Dental Laboratories

EMAIL

PHONE: **888.625.5757**

Know us by our Smiles!

[HOME](#) [PRODUCTS](#) [FORMS / SALES](#) [NEWS](#) [CAREER](#) [CONTACT](#)

Products

All Ceramic
Zirconia
PFM
Partials
Dentures
Implants

PFZ - Porcelain Fused Zirconia

Supreme accuracy in fitting in crowns and bridges – PFZ allows you to achieve an extraordinarily high degree of precision and accuracy in fitting.

Crown and bridge frameworks with a minimum wall thickness (0.4 mm) and interdental linkage (9 mm2) – as with high-strength dental alloys in the area of precious metals.



FCZ - Full Contour Zirconia

Full Contour Zirconia crowns are solid zirconia with no porcelain overlay. Your patients will be happy with a tooth shaded crown instead of metal occlusal PFM and full-cast gold dental restorations. Virtually chip-proof, FCZ solid zirconia crowns and bridges are the ideal dental restoration for your bruxing and grinding patients.



Zirconia Restored Implant

Please call for more information - **888.625.5757**

Exhibit 98

Tosch Laboratory, Inc.

A Certified CEREC inLab Crown and Bridge Dental Laboratory

LUNA Full Contour Translucent All Zirconia Crowns.


Luna translucent zirconia crowns are an esthetically pleasing alternative to full cast crowns, especially where strength is a concern. The Luna Zirconia crowns are virtually indestructible and therefore an excellent choice where there is minimal occlusal clearance or excessive interproximal space. Crowns and bridges are shaded prior to sintering so the color is internal and permanent.

"A recent study found that the degree of wear on the opposing tooth was actually four times greater with dental porcelain than in polished, unglazed zirconia, suggesting that zirconia restorations are not only more fracture-resistant than porcelain restorations, but also gentler in their functional interaction with natural dentition." In another study, "glazed Solid Zirconia crowns were found to wear compatible with enamel and virtually identical to glazed IPS e.max." Study copies are available at Tosch Laboratory.

Due to the strength of Luna Zirconia it is also an excellent choice instead of the traditional PFM bridge. It is possible to construct bridges as large as five and six units.

Contact us today to try a Luna Translucent Zirconia Crown

[Link to a study on the wear of full zirconia crowns versus standard dental porcelain on opposing dentition.](#)



[Return to Main Page](#)



Exhibit 99

THE ART OF SMILE TECHNOLOGY
CONTACT US **604.270.2230**


Solid Zirconia

Solid Zirconia Crowns and Bridges

These are suitable for both the anterior and posterior region. Buccal/Labial Cut-Back technique provides an excellent aesthetic appearance while maintaining high strength. Our zirconia crowns and bridges are designed and milled using CAD/CAM technology with the latest scanner from Sirona, inEOS Blue for consistent standardized quality.

Benefits:

- Biocompatible
- High strength
- No metal sensitivity
- Esthetical
- Low thermal conductible




— Solid zirconia Crown (before & after)

Indications:

- Bruxers and grinders who do not desire metal occlusal PFMs or full cast crowns
- Patients who desire a tooth colored restoration but have a lack of preparation space
- Patients with metal allergies

Preparation Requirements:

- Circular shoulder preparation with rounded inner edges or chamfer with a width of approx. 1 mm
- No feather edges
- Minimal occlusal reduction of 1.0mm
- Avoid sharp edges and angles
- Make sure to observe minimum layer thicknesses (stability of the restoration)




— Solid Zirconia crown with labial layering

Cementation Recommendation:

- Resin reinforced glass ionomer cement (RelyX Luting Cement, 3M ESPE; GC Fuji Plus, GC America)
- Resin cements for short or over-tapered preparations (RelyX Unicem, 3M ESPE; Panavia F2.0, Kuraray)

Exhibit 100



» Scanners

» Milling Machine

» Materials

ORIGIN Zirconia

Metals

Titanium

Cobalt-Chrome

Ceramics

Vitablocs®

e.max

Triceram®

Wax/PMMA/Temps

» Abutments

» CAM Software

» Sintering Oven


» Suction

» Consumables



- Strongest in Class (Documented 1,400 MPa)
- Tested clinically safe for cytotoxicity
- Far more chip resistant while finishing the sintered copings
- Excellent translucency, aesthetics and 100% biocompatibility
- High degree of clinical safety
- Engineered with bar codes for full control of consistent sintering shrinkage
- Conventional cementation
- Excellent marginal accuracy
- European CE & U.S. IdentCeram® certified

Premium Quality ORIGIN YZ Zirconia



Zirconia Advantages at a Glance

Zirconia is highly aesthetic

Zirconia is highly translucent especially when compared to traditional metal substructures and can be customized according to the 16 classic Vita shades. Zirconia restorations also solve the 'black gum line' problem inherent with traditional metal restorations which occurs when the gingiva recedes. These qualities of zirconia provide natural looking restorations which are indistinguishable from real teeth.

Zirconia is exceptionally strong

Often referred to as 'white steel', zirconia is just as strong if not stronger than metal restorations and has a high fracture resistance.

Zirconia is bio-compatible

The body accepts tissue friendly zirconia as a natural material which minimizes the potential for allergic reactions.

Zirconia is corrosion free

Zirconia does not corrode because it is "chemically inert".

Zirconia has low thermal conductivity


Most metals have a high thermal conductivity causing patients to feel greater sensitivity to fluctuations of hot or cold in the mouth. Zirconia restorations have a low thermal conductivity which helps minimize this problem.

Full Contour Solid Zirconia

The ORIGIN® FC Zirconia is a full contoured solid zirconia crown or bridge restoration with no porcelain overlay. ORIGIN FC Zirconia restorations are both a strong and highly aesthetic alternative to metal occlusal PFM and full-cast metal restorations. Virtually unbreakable and chip-proof, ORIGIN FC Zirconia crowns and bridges are the ideal restoration for posterior molar crowns when the patient desires a tooth colored restoration but lacks the preparation space for, or has broken a PFM crown in the past.

ORIGIN Full Contour Zirconia Advantages

- A more esthetic alternative to traditional full-cast or metal occlusal alloy restorations and is ideal for bruxers & grinders who have destroyed other restorations
- At 1132 MPa Origin® Full Contour Zirconia is strong enough for both single and bridge cases
- Made with biocompatible 100% Tosoh zirconia powder
- Conventional cementation using non-expanding resin modified glass ionomer
- Virtually unbreakable and chip-proof
- Fabricated using CAD/CAM technology for consistent fit and marginal integrity
- Backed by a 10 year free replacement guarantee




Send Us An Email
800.255.2839


Document Downloads

- ORIGIN CAD/CAM
- Intelligence Scanner
- 3Shape Scanner
- DuoTron Pro

[Back to Top](#)


[Back to Top](#)

The Strongest Zirconia



1. Source: DIFK Deutsches Institut für Feuerfest und Keramik GmbH, Bonn, Germany Material Strength of Zirconia Produced with Two Methods. Determined according to DIN EN 993-6 and quality assurance methods TOSOH. Date of testing: 8/19/2008
2. Source: www.wieland-dental-systems.com
3. Source: nobelbiocare.com
4. Source: www.stikon.de
5. Source: www.3m.com
6. Source: www.vita-zahnfabrik.com
7. Source: cercon-smart-ceramics.com

On the assumption of a maximum chewing strength of 290 N, it is postulated that an endurance limit of approximately 500 N is required, which is given at an initial strength of at least 1000 N. The characteristics of ORIGIN YZ zirconia can be summed up as follows:

ORIGIN YZ Zirconia	KIC=7.0MPa·m ^{1/2}	Fracture Toughness
	m=14;2	Weibull Parameter

Since the mid-1980s, zirconia has held a steady position in high-performance ceramics. More than 300,000 artificial hip joint implants and several thousand dental implants are impressive testimonies to the biological compatibility of zirconium dioxide.

There are no uncontrollable risks or long-term deficiencies to be expected in the area of dentistry associated with the use of high quality yttrium stabilized zirconia blocks.

- ORIGIN YZ Zirconia Discs Have Been Tested & Proven to Be Superior for Various Milling Tools and Methods with No Margin Chipping
- ORIGIN Has the Highest Fracture Toughness (KIC) of 7.0

ORIGIN - YZ Zirconia Disks

We have a solid and proven expertise regarding the relationship between blank density, the milling enlargement factor and the sintering shrinkage rate that lead to the consistency of the prosthesis from single copings to large bridge framework.

From this high-quality, yttrium stabilized zirconium oxide material you can make bridge frames up to 16 units or customized implant abutments.

ORIGIN YZ zirconia is biocompatible and particularly resistant to ageing. This material differs from conventional zirconium oxides in that it is extremely homogeneous. This is achieved by means of an advanced manufacturing process and results in an excellent fit, even for large bridgework.

Once you've chosen ORIGIN, you'll benefit from high quality zirconia at the most affordable material prices in the industry.

[Back to Top](#)

ORIGIN - YZ Zirconia Cube Blocks


Designed for use with Sirona's and MC XL machines, ORIGIN YZ zirconia milling blocks are used as a milled substructure using CAD/CAM technology both for veneering and the press-on technique. ORIGIN YZ zirconia is formulated using a small grain size (< 0.3µm) to reduce porosity and undergoes an advanced manufacturing process to ensure uniform density, no warping during the sintering process, and quality control for large ceramic blanks. The yttrium-stabilized high quality zirconium oxide blocks display exceptional physical properties including a flexural strength of 1100 MPa.

[Back to Top](#)

Tested and Proven to Be Superior for Various Milling Tools and Methods with No Margin Chipping


Margin chipping can occur for many reasons including; incorrect calibration, tool wear, a margin design that is too weak, etc. In addition, weak zirconia material may not withstand the stress of the aggressive milling process, possibly leaving small visible chips along the margin as well as invisible micro-cracks within the coping or framework.

ORIGIN YZ zirconia discs have been tested successfully for different milling tools including both conventional designs (tool B & C) and advanced burs (tool A) designed to cut faster and more aggressively. It also can be milled under both wet milling and dry milling conditions.

[Back to Top](#)

[Send Us An Email](#)
[800.255.2839](#)

Document Downloads

- [ORIGIN CAD/CAM](#)
- [Intelligence Scanner](#)
- [3Shape Scanner](#)
- [DuoTron Pro](#)



Excellent Esthetic Properties - A Translucent Substructure That is Comparable to or Better Than Leading Brands

An identical coping was milled with the Cerec machine using three of the major brands. The translucency of the material depends not only on the material properties of the ceramic, but also on the recommended thickness of the layer, i.e. the wall thickness. So we followed the minimum wall thickness requirement (0.7mm) of the inLab system for all three copings with the same design.

[Back to Top](#)

ORIGIN - YZ Zirconia Blocks and Disks are Engineered to be Used for the Following Milling Systems

- Cerec® inLab®
- Wieland Zeno®Tec
- KaVo Everest®
- AmannGirrbach Ceramill

[Back to Top](#)

16 Shades Possible - All Classic V Shades


Supports Implant Bar Design

- ORIGIN zirconia can be colored by dipping the coping into the Color Liquids prior to sintering.
- ORIGIN YZ zirconia blocks can be colored according to the 16 different colors covering the entire classic-V shade range.
- The color becomes indelible through the sintering process thus avoiding the possibility of having white margins.
- After dipping, the excess liquid is removed by blow drying.
- Coloring liquid has no negative effects on the physical material properties such as flexural strength or fracture toughness



[Back to Top](#)

[HOME](#) | [NEWS](#) | [PRODUCTS](#) | [DIGITAL SOLUTIONS](#) | [COMPANY](#) | [EDUCATION](#) | [EVENTS](#)
©Copyright B&D Dental Technologies



[Send Us An Email](#)
[800.255.2839](#)

Document Downloads

- [ORIGIN CAD/CAM](#)
- [Intelligence Scanner](#)
- [3Shape Scanner](#)
- [DuoTron Pro](#)

Exhibit 101

[Cart/Checkout](#) | [Reorder](#) | [Login](#) | [Register](#)
Search **3M ESPE**United States [\[Change\]](#)**Dental Products**

Buy it Now

[3M ESPE Dental Professionals](#) > [Products](#) > [Product Categories](#) > [Digital Materials and Equipment](#) > [Lava™ Zirconia](#)
Dental Products**Products**

- + Product Catalog
- Techniques and Solutions
- Brands You Trust
- Where to Buy 3M ESPE Dental Products
- Find an Authorized Lab

Promotions

- Special Offers
- Sign Up for 3M ESPE E-mail
- Digital Dentistry Newsletter

What's New

- New Dental Products
- News Articles
- Press Releases
- Upcoming Events

Resources

- Literature
- MSDS/Product Safety

Education

- Online Dental CE Courses
- Scientific Resources

Contact Us

- About Us
- Contact 3M ESPE
- Equipment Service


**Lava™ Zirconia**

- [Print](#)
- [Share](#)

[Overview](#) | [Features & Benefits](#) | [Indications](#) | [Technical Information](#) | [Ordering Information](#)

All ceramic crowns and bridges are an excellent long term esthetic solution for patients.

- High stability and precise fit
- Excellent long term natural aesthetics
- Biocompatibility



PFM

Zirconia

Like [Follow](#)**Customer Care**

1 (800) 634-2249

E-mail News

Stay in the know with exclusive e-mails from 3M ESPE. [Sign up »](#)

more information

Handling and Prep guide for Lava™ Restorations (pdf, 1.1MB)

Lava™ Zirconia – 10 Year Clinical Study (pdf, 1.8mb)

Why Lava™ Zirconia? (PDF, 957KB)

Lava™ Zirconia Scientific Facts (pdf, 1.3mb)

Lava™ Zirconia Scientific Facts - Zirconia-supported Ceramic Restorations: Uncovering the Mysteries (PDF, 765KB)

Lava™ Zirconia Scientific Facts - Zirconia is Not Alike (pdf, 354kb)

Lava™ Zirconia for Implant Abutments Brochure (PDF, 811KB)

Lava™ Zirconia for Implant Abutments Dentist Statement Stuffer (PDF, 512KB)

Lava™ Zirconia for Implant Abutments Dentist 30-second Sell Sheet (PDF, 564KB)

Lava™ Zirconia for Implant Abutments Handling Guide (pdf, 964kb)

EXHIBIT 101

-123-

Lava™ Zirconia Scientific
Facts – Implant Abutments
(pdf, 549kb)

Lava™
Precision Solutions

Read more about
Lava Materials and Equipment »

3M ESPE: [About Us](#) | [Site Map](#) | [Contact 3M ESPE](#) | [Equipment Service](#)

[Legal Information](#) | [Privacy Policy](#)

© 3M 2012. All Rights Reserved. **3M**



[Cart/Checkout](#) | [Reorder](#) | [Login](#) | [Register](#)
Search **3M ESPE**United States [\[Change\]](#)**Dental Products**

Buy it Now

[3M ESPE Dental Professionals](#) > [Products](#) > [Product Categories](#) > [Digital Materials and Equipment](#) > [Lava™ Zirconia](#)
Dental Products**Products**

- + Product Catalog
- Techniques and Solutions
- Brands You Trust
- Where to Buy 3M ESPE Dental Products
- Find an Authorized Lab

Promotions

- Special Offers
- Sign Up for 3M ESPE E-mail
- Digital Dentistry Newsletter

What's New

- New Dental Products
- News Articles
- Press Releases
- Upcoming Events

Resources

- Literature
- MSDS/Product Safety

Education

- Online Dental CE Courses
- Scientific Resources

Contact Us

- About Us
- Contact 3M ESPE
- Equipment Service

**Lava™ Zirconia**

- [Print](#)
- [Share](#)

[Overview](#) [Features & Benefits](#) [Indications](#) [Technical Information](#) [Ordering Information](#)
All Zirconia is not alike!

Not all Zirconia is alike. Yes, Zirconia can be chemically similar, but once processed the mechanical and optical properties can be significantly different. Lava™ Zirconia has now been commercially available for 10 years.

• **Productivity**


- High strength
- High stability

• **Excellent long term aesthetics**

- Tooth colored material
- Unique shading system
- Homogeneous coloring
- High, life-like translucency

• **Biocompatibility**

- Metal free
- Clinical studies showed no allergic reactions (gingival)
- Excellent regeneration of gingiva

• **Excellent fit**Like [Follow](#)**Customer Care**

1 (800) 634-2249

E-mail News

Stay in the know with exclusive e-mails from 3M ESPE. [Sign up »](#)

more information

Handling and Prep guide for Lava™ Restorations (pdf, 1.1MB)

Lava™ Zirconia – 10 Year Clinical Study (pdf, 1.8mb)

Why Lava™ Zirconia? (PDF, 957KB)

Lava™ Zirconia Scientific Facts (pdf, 1.3mb)

Lava™ Zirconia Scientific Facts - Zirconia-supported Ceramic Restorations: Uncovering the Mysteries (PDF, 765KB)

Lava™ Zirconia Scientific Facts - Zirconia is Not Alike (pdf, 354kb)

Lava™ Zirconia for Implant Abutments Brochure (PDF, 811KB)

Lava™ Zirconia for Implant Abutments Dentist Statement Stuffer (PDF, 512KB)

Lava™ Zirconia for Implant Abutments Dentist 30-second Sell Sheet (PDF, 564KB)

Lava™ Zirconia for Implant Abutments Handling Guide (pdf, 964kb)

EXHIBIT 101 -125-

Lava™ Zirconia Scientific
Facts – Implant Abutments
(pdf, 549kb)

Lava™
Precision Solutions

Read more about
Lava Materials and Equipment »

3M ESPE: [About Us](#) | [Site Map](#) | [Contact 3M ESPE](#) | [Equipment Service](#)

[Legal Information](#) | [Privacy Policy](#)

© 3M 2012. All Rights Reserved. **3M**



[Cart/Checkout](#) | [Reorder](#) | [Login](#) | [Register](#)Search **3M ESPE**United States [\[Change\]](#)**Dental Products**[Buy it Now](#)[3M ESPE Dental Professionals](#) > [Products](#) > [Product Categories](#) > [Digital Materials and Equipment](#) > [Lava™ DVS Full Contour Digital Crowns](#)**Dental Products****Products**

- Product Catalog
- Techniques and Solutions
- Brands You Trust
- Where to Buy 3M ESPE Dental Products
- Find an Authorized Lab

Promotions

- Special Offers
- Sign Up for 3M ESPE E-mail Digital Dentistry Newsletter

What's New

- New Dental Products
- News Articles
- Press Releases
- Upcoming Events

Resources

- Literature
- MSDS/Product Safety

Education

- Online Dental CE Courses
- Scientific Resources

Contact Us

- About Us
- Contact 3M ESPE
- Equipment Service

Lava™ DVS Full Contour Digital Crown• [Print](#)• [Share](#)

Overview	Clinical Cases	How to Order
----------	----------------	--------------

Lava™ DVS full contour digital crowns deliver the strength, fit and beauty you and your patients require.

Now your patients can experience the beauty of a digital Lava restoration that is more esthetic than pressed ceramics. And the digital precision of a Lava DVS full contour crown means that manual adjustments are minimized and fit is optimized.



Final restoration after characterization

More esthetic than pressed ceramics

Left: Lava DVS Zirconia Crown; Right: Pressed Zirconia Crown

- Available in a variety of **shades and translucencies**. You're sure to receive the restoration that closely mirrors your patient's natural esthetics.
- Indicated for both **posterior and anterior** restorations.
- For use with **conservative prep** and conventional cementation.
- Digital precision** ensures consistent strength, minimal adjustments and optimized fit.

Consistent Strength

The strength of Lava™ DVS crowns makes them an ideal solution. Why? Our controlled CAD/CAM process and a strong zirconia core ensure Lava DVS restorations are strong and beautiful every time.

Like [Follow](#)**Customer Care**

1 (800) 634-2249

E-mail News


Stay in the know with exclusive e-mails from 3M ESPE. [Sign up »](#)

Lava™**Precision Solutions**

Read more about
[Lava Materials and Equipment »](#)

More Information

Lava™ DVS Digital Veneering System Brochure for the Dentist (PDF, 563KB)




3M ESPE: [About Us](#) | [Site Map](#) | [Contact 3M ESPE](#) | [Equipment Service](#)

[Legal Information](#) | [Privacy Policy](#)

© 3M 2012. All Rights Reserved. **3M**



Exhibit 102



Dental Prosthetic Services

A Leader in Quality and Customer Service Since 1972


[Home](#)
[About Us](#)
[Products](#)
[Dental Sleep
Medicine](#)
[Services](#)
[Events](#)
[Downloads](#)
[Employment](#)
[Contact Us](#)

Downloads

[Doctor Preference Form](#)

[Smile Guide](#)

[RX forms](#)

[Request Rx Pads and Boxes](#)

[Join Our Mailing List](#)



Testimonials

"Whether it is implant dentistry, cosmetic veneers, dental appliances, or crown and bridge cases, DPS has consistently provided our dental team and patients with quality service and experience for many years."

- Dr. Jon R. Hart
- Dr. Peter J. Blodgett
- Dr. John E. Wilmes

"DPS plays an integral role in the exceptional care we provide for our patients."

- Dr. Patricia Meredith

Like all zirconia restorations, our Nautilus restorations offer you superior strength. But, with its enhanced anatomy and availability in Classic or 3D shades, our Nautilus Full-Contour offers you a more esthetic all-zirconia option in the posterior. And, our Nautilus Core combines the strength of a zirconia coping with the beauty of a stacked porcelain, rivaling all-ceramic and PFM alternatives.

Nautilus Full-Contour

Nautilus Full-Contour is now better than ever!

All Nautilus Full-Contour restorations designed after January 18 will have our new enhanced anatomy.



Nautilus Full-Contour crowns and bridges are solid zirconia restorations with no porcelain overlay. At 850 MPa, Nautilus Full-Contour is one of the strongest restorations available. Our technicians design Nautilus Full-Contour with enhanced anatomy and shade it to match your patient's natural tooth color; giving you a more esthetic all-zirconia option.


Ideal for posterior crowns and bridges (up to 8-unit spans), patients now have a strong alternative to full-cast restorations. Nautilus Full-Contour can also be used in the anterior for patients who need a stronger anterior restoration.

And while it is incredibly strong, polished zirconia is significantly less abrasive than standard porcelain, making it kinder to the opposing dentition than other restorative ceramic material.

Prep Requirements

- A chamfer margin is preferred; however any traditional full-cast crown prep is acceptable. Avoid a beveled or knife-edge design.

Exhibit 103



Call: 905 665 9990 1 855 MY 5 AXIS

1380 Hopkins St. Unit 1 Whitby ON Canada L1N 2C3

5 Axis Dental Design Center

[Home](#) [Sending Cases](#) [About Us](#) [Print Forms](#) [iTero](#) [Products](#) [Our Technology](#) [Links](#) [Contact Us](#) [Materials](#)



Full Zirconia Crowns


57573

PATENT

Zirconia full anatomical crowns are milled using several different sizes of precision burs that allow us to create fine occlusal details. Recommended for heavy bruxers for posterior restorations.

Manufactured from Metoxit TZP-A zircona ($ZrO_2/Y_2O_3/Al_2O_3$). With a compressive strength of 2000 MPa and a flexural strength of 1200 MPa. Metoxit AG Switzerland brand zirconia - FDA and Health Canada approved.





Developed By: **THE ONLINE PRACTICE**

Whitby ON Canada Dentist | Full Zirconia Crowns. Dr. Pamela Mursic is a Whitby ON Canada Dentist.

5 Axis Dental Design Center is not affiliated with any of the following companies:

Lifecore Biomedical, Inc., Biomet 3i, Nobel Biocare

Steri-Oss, Stryker, Sustain, APM Stemgold, Preat, Soodenco, Parkell, O-IC Implant, IMTEC Corp, Smooth Staple, Stemgold-ImplaMed, Cendres & Metaux, Implant

Direct, etc. Ankylos® is a federally registered trademark of DeguDent GmbH, Germany

Astra® is a federally registered trademark of Astra Tech AB, Sweden

Branemark®, Steri-Oss®, Steri-Oss Replace®, Replace Select®, Steri-Oss Hex-Loc®, NobelReplace®, NobelSpeedy®, NobelPerfect®, NobelDirect® are federally

registered trademarks of Nobel Biocare AG, Zurich, Switzerland. BioHorizons® is a federally registered trademark of BioHorizons Implant Systems, Inc. Mini-SC®, Dolder®, Biloc®, Combi-Snap®, Dalbo® are federally registered trademarks

of Cendres & Metaux, Switzerland

GI-MASK® is a federally registered trademark of Coltene/Whaledent AG, Switzerland

Entegra®, Endopore®, Pitt-Easy® and Anatomic® are federally registered trademarks of Innova LifeSciences Corp. Toronto, ON Canada 3i®, Osseotite®, Osseotite Certain®, Miniplant® and MicroMiniplant® are federally registered

trademarks of Biomet 3i, Inc. Florida, USA IMZ® is a federally registered trademark of IMZ GmbH, Germany

ITI® Dental Implant System, Straumann® and synOcta® are federally registered trademarks of Institut Straumann, Switzerland.

Legacy® is a federally registered trademark of Implant Direct LLC, USA. Lifecore®, Restore®, Sustain®, Stage-1®, Renova® and Prima® are federally registered trademarks of Lifecore Biomedical, Inc., Chaska, MN. Xive TG® is a federally registered trademark of Friadent GmbH, Germany. Zimmer Dental®, Screw-Vent®, Core-Vent®, Bio-Vent®, Micro-Vent®, Spline®, Advent®, SwissPlus®, Integral®, Omniloc®, Threadloc®, Zimmer One-Piece® are federally registered trademarks of Zimmer Dental, Inc., Carlsbad, CA.

Exhibit 104



(877) 954-6243

Virtual Smiles

My Account

[HOME](#) | [PRODUCTS](#) | [GALLERY](#) | [TESTIMONIALS](#) | [ABOUT US](#) | [CONTACT US](#)

- [Zirconia Crowns](#) ›
- [IPS e.max](#) ›
- [PFM](#) ›
- [Implants](#) ›
 - [Custom Implant Abutments](#)
 - [Screw-retained Denture](#)
 - [Hybrid Implant Bridge](#)
 - [Overdentures](#)
- [Milling Center](#) ›
- [Removables](#) ›
 - [Bayflex](#)
 - [Complete Dentures](#)
 - [Partials](#)
- [Fierce Mouthguards](#) ›
- [Orthodontics](#) ›
 - [Retainers & Finishing](#)
 - [Study Models](#)
 - [Space Maintainers & Regainers](#)
 - [Arch Development](#)
 - [Functional Appliances](#)
 - [Molar Distalizing](#)
 - [Splints](#)
 - [Holding Appliances](#)
- [Other Services](#) ›

Zirconia Crown Dental Lab


Top Quality Zirconia Crowns

Bayshore offers a variety of crown options using only the most technologically advanced, highest quality materials available to fit your needs. Our highly trained crown dental lab technicians will guide you in choosing the best option for each individual case, no matter the complexity of the restoration.

Some of our most popular products are the Z Crown™ Solid Zirconia Crown, and the Z Crown™ Layered Zirconia Crown.

Z Crown™ Solid

Full Contour Zirconia Crown



Our Z Crown™ Solid is made from yttria stabilized Zirconia, creating an esthetic, biocompatible, and completely chip-resistant crown with a natural, vibrant translucency. Fabricated using Bayshore's state-of-the-art CAD/CAM technology, quality and consistency are ensured in every case.

**Make the move to metal free
without compromising strength,
quality, or consistency.**

The Z Crown™ Solid is *the* esthetic alternative to full cast crowns.

Features

- On-site, state-of-the-art CAD/CAM technology
- 1200 MPA flexural strength
- High level of biocompatibility

Benefits

- Esthetic alternative to full cast crowns and metal occlusals
- Completely chip-resistant (solid Zirconia)
- Saves tooth structure (minimal shoulder preparation)
- Glazed final restoration prevents plaque buildup and opposing dentition wear
- Precise marginal fit

If you would like more information on Z Crown™ Solid, or our full-service crown dental lab, [contact us](#) today!

Exhibit 105




All Products Made in the US

Calls us at (800) 243-2235
(781) 438-6110

[HOME](#)
[PRODUCTS](#)
[DIGITAL TECHNOLOGIES](#)
[CONTACT](#)
[RESOURCES](#)
[BLOG](#)
[EVENTS](#)
[ABOUT US](#)

ADT Full Zirconia

High Translucency Full contour Zirconia Restorations




ADT FZ is milled from the most esthetic full contour zirconia available. The Material is isostatically pressed zirconia blocks manufactured in Germany. We mill your monolithic restoration and perform a multi-color shading process that dramatically enhances the esthetics.

ADT FZ is indicated for full crowns and bridges up to 12 units. This product is available in 16 Vita shades. To enhance the esthetics even further, we offer five modifier shades: Violet, Gray, Brown, Orange and Pink. Tooth preparation requirements can be that of a PFM or an all ceramic material. In addition, full contour zirconia can be conventionally cemented.

In addition to its strength and beauty, multiple studies have shown full zirconia restorations to be gentler to surrounding dentition than conventional PFM's.

- Lifetime guarantee against breakage
- Chip Resistant
- Less wear to adjacent teeth than PFM
- Conventional cement protocol
- Conventional tooth preparation
- 1100 Mpa flexural strength
- Three times the strength of e.max
- Available in 16 Vita shades with
- Competitive pricing helps increase your profitability
- Award winning material (JDT 2011 WOW Product of the year!)




ADT's technical capabilities ensure these products are milled precisely to provide you with excellent fitting restoration for fast seating times. The combination of automation and use of a low cost, yet high quality material allows us to offer you the best restorations for your patients while saving you time and money.

Indications

Full contour zirconia for crowns and bridges up to 12 units. Ideal product for replacement of full gold crowns and PFMs.

Technical Specs

- Flexural strength - 1,100 MPa
- Compressive Strength - 3,000 MPa
- Conventional cementation
- Conventional tooth preparation



Cementation

- RelyX™ Ultimate Adhesive Resin Cement with Scotchbond™ Universal Adhesive
 - RelyX™ Unicem 2 Self-Adhesive Resin Cement
- Please see the manufacturer's Instructions for Use for the proper cementation procedure.

ADA Codes

D2740 Crown

D6245 Pontic Porcelain/Ceramic

Links to wear studies of FZ to enamel:

A study on the in-vitro wear of the natural tooth structure by opposing zirconia. Journal of Advanced Prosthodontics.

Wear Study of Zirconia to opposing enamel. Clinicians Report.

Two-body wear of different ceramic materials opposed to zirconia ceramic. Journal of prosthetic Dentistry

Ask Dr. Christensen


Our Full contour Zirconia vs the sledge hammer!

As an example of the strength of our full contour zirconia we tested it (and other materials) against a 14 pond sledge hammer. For results click on the arrow to the right.



© 2012 Advanced Dental Technologies

Exhibit 106



800.337.7800



PRODUCTS

DENTIST SUPPORT

SERVICE CENTER

TESTIMONIALS

ABOUT US

CONTACT US

PH: 877.337.7800



Custom-crafted to your specifications.

(And we'll sink our teeth into any challenge.)

PFM Full Cast Restoration

All Ceramics

Cercon

IPS Empress

Wol-Ceram

IPS e.max

Z-Crown Solid Zirconia

Z-Crown Layered Zirconia

Removables

Implants

Fierce Mouthguards

Orthodontics

Other Services & Supplies

Z-Crown Solid Zirconia

Z-Crown Solid is monolithic and made from 100% pure zirconia. At 1200 MPa, the Z-Crown has a comparable flexural strength of traditional porcelain-fused-to-metal restorations while maintaining a natural and vibrant translucency.

The product 100% metal-free, a feature that prevents gingival darkening and removes the possibility of exposing metallic margins if gum recession begins. They are naturally esthetic, transmit the color of adjacent teeth and can be matched to any shade. Z-Crown Solid is one of the strongest restorations available for any patient who may have bruxing or grinding issues. Utilizing CAD/CAM technology, it also has a precise marginal fit that can significantly reduce chair time. Since traditional cementation is recommended, there is also no need to bond the restoration.

Features

- 1200 MPA flexural strength
- Highly biocompatible
- Low CTE
- CAD/CAM technology
- CAM only option
- White and shaded blanks
- Any design
- Non-Feldspathic porcelains
- Bond or cement

Benefits

- Long-term reliability
- Healthy tissue response
- Synthetic ceramics for low wear
- Laser precision
- Custom position pontics — without library
- Better shading
- Same as alloy without special needs
- Natural opalescence and translucency
- Versatility of choice

At DDS Lab, we have technicians on staff whose expertise aided the research and introduction of the Z-Crown product to the dental profession in the United States. Our technicians can create crowns and bridges of any size and specification with laser precision and esthetic luminescence.



Z-Crown Zirconia Product Gallery



I'm shipping my case on:

Date:



"The technicians at DDS Lab truly have become a part of my team. They follow my prescription and diligently work with me on my cases..."

*especially the really
complicated cases."*

*-Dr. Susannah Bohannon
Olympia, WA*



©2011 DDS Lab. All Rights Reserved | [Contact Us](#) | [Legal](#) | [Privacy Policy](#)


Exhibit 107

[HOME](#) [LOGIN](#) [REGISTRATION](#)


[• Articles/Studies](#) [• News](#) [• Regulatory](#) [• Testimonials/Videos](#)

***"Because every child
is special!"***

...And So is Every Smile."


[Product List](#) | [Order Online](#) | [Contact Us](#) | [Related Pediatric Sites](#)


NuSmile ZR represents the perfect balance between art and science. Made from high grade monolithic zirconia ceramic, NuSmile ZR looks like the real thing and is 9 times stronger than natural teeth. Easy to place, NuSmile ZR offers the highest esthetics to satisfy the most demanding clinicians and parents. Our latest crown offers superior durability to last until the tooth exfoliates without needing any ongoing repairs. The translucency of our zirconia ceramic has been optimized for both natural esthetics and to prevent the problem of dark tooth show through of pulpally treated teeth, ensuring continued satisfaction with the restoration over time. **Nothing beats a NuSmile.**

[General Information](#)
[Case Examples](#)
[Instructions for Use](#)
[Watch ZR](#)
[Video/Request Sample](#)


Superior esthetics
Anatomically contoured
Strong and metal free
Easy to use - conservative preparations
Superior durability
Excellent biocompatibility
Autoclavable

Follow this link to our sister company:



laboratory specializing in orthodontic,
pedodontic and TMD appliances.

1-800-346-5133



©2004-2012 NuSmile® Pediatric Crowns. No further publication or commercial use may be made of the materials on this website without the express written consent of NuSmile® Pediatric Crowns. All rights reserved.

Exhibit 108




Welcome to Crown Dental
Laboratory!


(800)680-8608

Home
About Us
Dentures and Partial
Crown and Bridge
Implants
Prettau Bridge
Staub Cranial Technique
About our Leadership
Lab Schedule
Contact Us
Home
About Us
Dentures and Partial
Crown and Bridge
Implants
Prettau Bridge
Staub Cranial Technique
About our Leadership
Lab Schedule
Contact Us

Crown Dental Lab can meet
all of your implant
restoration needs, from
Removable solutions like
Locator Over-dentures, to
the absolute latest in Screw
Retained Hybrid
Restorations like our full
zirconia Prettau Bridge.




Implant Locator Over-Dentures



360Zr Full Contour Zirconia
Crown

Let Crown Dental Lab show
you how you can provide a
strong, full contour
Zirconia Crown at a price
that will meet every budget.
Ask us about the 360Zr
Crown.





Prettau Bridge Full Zirconia
Hybrid

Crown Dental Lab hand mills
this full arch out of one solid
disc of Zirconia. This screw
retained Full arch Hybrid
gives you the strength and
esthetics that will excite your
patients. Learn more about
the Prettau Bridge from
Crown Dental Lab.

STAUB
CRANIAL

**Ask us about our
expertise in** Screw
Retained Hybrid
Restorations



*Deliver dentures in just 3
appointments!!* Let Crown
Dental Lab, a fully certified
Staub Cranial lab, show you
how easy it can be.

**Check out our
turnaround times by
clicking above!**

Website powered by Network Solutions®




Home
About Us
Dentures and Partial
Crown and Bridge
Implants
Prettau Bridge
Staub Cranial Technique
About our Leadership
Lab Schedule
Contact Us
Home
About Us
Dentures and Partial
Crown and Bridge
Implants
Prettau Bridge
Staub Cranial Technique
About our Leadership
Lab Schedule
Contact Us


Crown and Bridge

Our Crown and Bridge department offers a full range of choices to meet the specific treatment plans you have chosen for your patients. We provide traditional PFM in High Noble, Semi-precious, Non-precious, Zirconia and Titanium.

If your patient prefers the option of full metal restoration, we can provide you with the metal solution of your choice. Or, if you prefer a more natural appearance, we offer a full contour, solid Zirconia option in any Vita Shade required.



360Zr Full Contour Zirconia
Crown




Full Cast Crowns



Multi-Unit Bridge

If you would like more details regarding the full capabilities of our Crown and Bridge Department, please contact us by email at: mikestafford@crowndental.biz or call (800)680-8608 and ask for Mike Stafford.



Anterior PFM Crowns



Diagnostic Wax Up

Website powered by Network Solutions®

Exhibit 109


[Products](#)[Back to the product gallery](#)[CAP FZ - Full Zirconia](#)

CAP FZ - Full Contour Zirconia

[Introduction](#)[Features and Benefits](#)[Stain and Glaze Instructions](#)[Resources and Downloads](#)[Technical Specs](#)[Request a Sample](#)[Lifetime Warranty](#)

CAP FZ - Full Contour Zirconia with MultiShade

CAP FZ is milled from the most esthetic full contour zirconia available. We mill your restoration and perform a multi-color shading process that dramatically enhances the esthetics. We then sinter the restoration and return it to you ready to stain and glaze. CAP FZ is indicated for full crowns and small bridges. We are now accepting digital files from 3Shape or Dental Wings, or we can produce CAP FZ from your models using our proprietary anatomic library. We will return your restoration ready for stain and glaze.



Journal of Dental Technology's
WOW!
2011

IDENTCERAM

Custom

CAP FZ |



08

**Lifetime
Against**See terms: <http://cap-us.com/materials/details/28>
Denti

Identceram stickers shipped
with every CAP Full Zirconia
order

CAP Warranty stickers
shipped with every qualifying
order

Learn why CAP FZ is better and how it's made



CAP Warranty

Pellentesque habitant morbi tristique senectus et netus et malesuada fames ac turpis egestas. Vestibulum tortor quam, feugiat vitae, ultricies eget, tempor sit amet, ante. Donec eu libero sit amet quam egestas semper. Aenean ultricies mi vitae est. Mauris placerat eleifend leo.

HOW TO BUY

FIND OUT MORE


Exhibit 110



TOMORROW'S TECHNOLOGY...
YESTERDAY'S PERSONAL SERVICE

HOME
PRODUCTS/SERVICES
CUSTOMER CENTER
CE COURSES
CONTACT US
ABOUT US

WORK
AUTHORIZATION
IVOCLEAN
BROCHURE



What is TLZ?

TLZ is our all-zirconia, monolithic crown or bridge restoration. TLZ is guaranteed against breakage for the life of the patient. Tested at a flexural strength of 850 MPa, our High-Translucency, monolithic formula is stronger than nearly all other restorations on the market. A unique fabrication process makes TLZ beautiful, lifelike and a high recommendation from ROE.




TLZ is appropriate for posterior restorations, is completed in just 6 laboratory days and is guaranteed to satisfy.

When you need the strongest restoration available, ask for TLZ!

Crowns for Life™

Limited Lifetime Warranty Statement

Crowns for Life™ applies to ROE Dental Laboratory's TLZ all zirconia crowns. TLZ is an acronym for Translithic Zirconia, the material used to fabricate TLZ crowns. This is a warranty against breakage resulting from normal mastication. This is not a warranty against shade change, wear, or replacement required from additional dental procedures. This limited warranty covers replacement of the TLZ crown only, and not any costs associated with removal or replacement. Replacement with an alternate restorative material may result in additional costs.



How is it designed for dental?

Zirconia, which is unearthed as zirconium dioxide, is a white crystalline oxide that exists in three phases, cubic, tetragonal, and monoclinic, depending on temperature and pressure formation. In dentistry we concern ourselves with the latter two. Zirconia in the pure tetragonal phase is unstable. In order to create the milling blocks used in our machines, dental manufacturers add yttria, creating yttria-stabilized tetragonal polycrystals (Y-TZP)₂. In this state zirconia is extremely hard and possesses a unique characteristic called transformation toughening. When tensile stress is introduced from crack propagation, the tetragonal formation morphs into monoclinic, increasing the volume 3 to 5%, and subsequently transforms the stress from tensile to compressive. This self-healing mechanism makes zirconia the ideal material for the oral environment.

In-operative adjustment of zirconia As mentioned above, micro crack propagation is prevented while the monoclinic transformation occurs, thereby increasing surface tension and tensile strength (i.e., - theoretically, grinding can

Cerec In-Lab / Cadent
iTero
Porcelain To Zirconia
Monolithic Zirconia
IPS Empress Esthetic
IPS E.max
TLZ Zirconia
Lava
Bellglass
Tescera
PFM & All Cast
Other Fixed Restorations

increase the strength of Y-TZP zirconia). But before grinding-away on your crowns, you should observe certain criteria. The severity of grinding and the rise in temperature will affect the volume percentage of toughening. We recommend using fine-grade burs and copious amounts of water coolant to decrease heat generation. Researchers have studied the quality of diamond particles impregnated in dental burs, the hardness of the binding material, and the precision and centricity of the shafts, and have found that no one manufacturer has a superior product. All manufacturers agree that hydration and very light-to-no pressure is the best technique – let the bur do the work and avoid dull tools. Also, the inside of a coping or monolithic restoration must be left untouched. It is recommended that if internal adjustments are needed for seating, the preparation should be adjusted.

Wear to opposing A common question that arises about our TLZ allzirconia restorations concerns the wear effects they have on opposing dentition. There are many documented studies measuring the wear of zirconia against fluorapatite, porcelain, gold, lithium disilicate, lucite etc. These invitro studies can be summarized by focusing on particle size and finished surfaces. Zirconia is comprised of ultra-fine particles that do not become saw-tooth when roughened, unlike standard porcelains used in dental restorations. This characteristic evidently leads to low-wear of opposing enamel. A Study in J Adv Prosthodont 2010;2:111-5 showed that wear to the antagonist teeth is much less than that of feldspathic porcelains. Moreover, the study agrees with many other studies we've reviewed, zirconia shows very low wear when highly polished.

Rella Christiansan's TRAC Research group has released preliminary results of a 7-year full-contour zirconia wear study which supports both of these claims. The study seeks to measure the amount of wear zirconia and other monolithic restorations exhibit invivo. After one year, the "very promising" results show that zirconia "mimics" natural dentition. Currently many other companies are completing invivo wear studies. We will keep you informed as further studies are published.

ROE RECOMMENDATION

Posterior crowns

PRIMERS & CLEANING AGENT:

Use primers when mechanical retention is minimal. They generate a sound bond between the restoration and the cementation material Monobond

Ivoclean™ & Monobond™ Plus - Ivoclar Vivadent
Monobond Plus - Ivoclar Vivadent
ZPrime - Bisco
Scotchbond - 3M

FLEXURAL STRENGTH

850 mpa

PREPARATION

Conventional C&B(chamfer)

RECOMMENDED BUR KITS FOR ADJUSTING ZIRCONIA:

SS White® The Great White® Z
Komet USA LD0540 Zirc / emax Adjustment Kit
Brassler USA KO196

ADA CODES

D2740 Crown - Porcelain/Ceramic Substrate
D6245 Pontic Porcelain/Ceramic
D6740 Abutment Crown Porcelain/CeramicC



Glidewell Laboratories v. Keating Dental Arts, Inc.
United States District Court, Central, Case No. SACV11-01309-DOC (ANx)

CERTIFICATE OF SERVICE

I hereby certify that on November 26, 2012, I electronically filed the document described as **SUPPLEMENTAL APPENDIX OF EVIDENCE IN SUPPORT OF GLIDEWELL'S OPPOSITIONS TO KEATING'S MOTIONS FOR SUMMARY JUDGMENT** with the Clerk of the Court using the CM/ECF System which will send notification of such filing to the following:

David G. Jankowski
Jeffrey L. Van Hoosear
Lynda J Zadra-Symes
Knobbe Martens Olson and Bear LLP
2040 Main Street, 14th Floor
Irvine, CA 92614

**Attorneys for Defendant Keating
Dental Arts, Inc.**
Tel: (949) 760-0404
Fax: (949) 760-9502

Jeffrey.vanhoosear@kmob.com
David.jankowski@kmob.com
Lynda.zadra-symes@kmob.com
litigation@kmob.com

Dated November 26, 2012

SNELL & WILMER L.L.P.

By: s/Philip J. Graves

Philip J. Graves
Greer N. Shaw

Attorneys for Plaintiff
James R. Glidewell Dental Ceramics, Inc.
dba GLIDEWELL LABORATORIES

16139994.1

Certificate of Service
SACV11-01309-DOC (ANx)