

Dental Zirconia Crown - BruxZir Solid Zirconia Crowns & Bridges - Glidewell Direct BruxZir Dental Zirconia Lab

 **GLIDEWELL LABORATORIES** [HOME](#)

DENTIST LABORATORY PATIENT **Product FOCUS** BRUXZIR IMPLANTS SNOORING & APNEA MY ACCOUNT

SERVICES | PRODUCTS | CAD/CAM SYSTEMS | EDUCATION | VIDEO GALLERY | LAB PERSPECTIVES | PARTNERSHIP PROGRAMS

DOWNLOAD RX

ALL-CERAMICS

Prismatik Clinical Zirconia

BruxZir Solid Zirconia

Lava Crowns & Bridges

Vivaneers No-Prep Veneers

IPS e.max & IPS Empress

DIAGNOSTIC WAX-UPS

IMPLANTS

PFMS

COMPOSITES & PROVISIONALS

REMOVABLES

BruxZir® Solid Zirconia

 Product  Technical  Video



BruxZir® Solid Zirconia is a full contour solid zirconia crown or bridge restoration with no porcelain overlay. More brawn than beauty, BruxZir is an esthetic alternative to posterior metal occlusal PFM and full-cast metal restorations. Virtually chip-proof, BruxZir Solid Zirconia crowns and bridges are the ideal restoration for bruxers and grinders.

If you would like to market BruxZir Solid Zirconia to your dentists, we offer three ways for you to become an Authorized BruxZir Lab.

AVAILABLE SERVICES

GLIDEWELL DIRECT
888-303-3975 U.S.
888-278-0414 Canada

Crown or Bridge on Model
Crown or Bridge from Scan
Mill Crown or Bridge on
CEREC InLab MCXL
Mill Crown or Bridge on Round Blank

PDF DOCUMENTS

BruxZir vs. Ceramco 3 - A Comparative Wear Study
BruxZir Enamel Wear Test Results
BruxZir Study: Promising 1-year Results
Instructions for Adjusting and Polishing BruxZir crowns & bridges
Do not use discs to finish full-contour zirconia

EXHIBIT *q3*
WIT: *Portelo*
DATE: *10-23-12*
L. Moskowitz, CSR 10816, RPR

Dental Zirconia Crown - BruxZir Solid Zirconia Crowns & Bridges - Glidewell Direct BruxZir Dental Zirconia Lab

If you own a 3Shape Dental System scanner with 2009 DentalDesigner™ software, upload your 3Shape design file to Glidewell Dental Lab. We mill — you color, sinter, stain and glaze. Or, pay \$10 more and we will use your 3shape design file to mill, sinter and color your BruxZir Solid Zirconia restoration — you stain and glaze.

You may also send your master & opposing modelwork with pinned, trimmed dies and a bite registration. Then choose between staining and glazing BruxZir restorations at your dental lab or to receive finished BruxZir Solid Zirconia restorations.

If you own a milling system, you can purchase the special BruxZir high translucent dental zirconia blanks. The 98 mm diameter BruxZir zirconia blanks are available in four thicknesses: 12 mm, 15 mm, 20 mm and 25 mm.

BruxZir is made from translucent PrismaTik Clinical Zirconia™. BruxZir is a registered trademark of Glidewell Laboratories. PrismaTik Clinical Zirconia is a trademark of Glidewell Laboratories.

Instructions for Seating BruxZir and
other zirconia-based crowns & bridges

ADDITIONAL RESOURCES

BruxZir Authorized Lab Partnership
Program

EXTERNAL LINKS

www.BruxZir.com

[BruxZir Blog](#)

[BruxZir Twitter](#)



[Shipping](#) | [Company](#) | [Careers](#) | [Contact Us](#) | [Conventions](#) | [Downloads](#) | [Site Map](#)

[Policies & Warranty](#)



©2012 Glidewell Laboratories
4141 MacArthur Blvd, Newport Beach, CA 92660
800-854-7256 USA | 888-278-0414 Canada



aveoTSD



BruxZir



Inclusive



FastScan



H/S Splint



Facebook



Twitter



iTunes



YouTube